

Surgical Blades and Scalpels

Surgical Blades and Scalpels

KYUAN®



Production Capacity can reach to 100 million pieces per year

No.10		No.20	
No.10A		No.21	
No.11		No.22	
No.12		No.23	
No.14		No.24	
No.15		No.25	
No.15C		No.36	
No.18			

■ Sterile Surgical Blades

- ▼ Available in both Stainless Steel and non-glare Carbon Steel.
- ▼ The best raw material and the production line of the latest technology ensure the quality and production capacity.
- ▼ Individually wrapped in sterile aluminum foil pack by Gamma Radiation.
- ▼ 100 blades in a dispenser box. 5000 pcs per carton.
- ▼ Carton Size(M): 0.37*0.20*0.16
- ▼ Gross Weight: 7.5KG/Carton
Net Weight: 7.0KG/Carton

■ Sterile Disposable Scalpels

- ▼ Available in both Stainless Steel and non-glare Carbon Steel.
- ▼ The blades are safely protected by LDPE guard.
- ▼ Individually wrapped in sterile aluminum foil pack by Gamma Radiation.



No.3 Scalpel


- ▼ 10 pcs in a dispenser box, 500 pcs/Carton.
- ▼ Carton Size(M): 0.43*0.30*0.20
- ▼ Gross Weight: 6.5KG/Carton Net Weight: 6.0KG/Carton



No.4 Scalpel

■ Stitch Cutters

- ▼ Used for general suture removal and designed with a rounded safety and bull nose tip.
- ▼ Available in both Stainless Steel and non-glare Carbon Steel.
- ▼ Individually wrapped in sterile aluminum foil pack by Gamma Radiation.
- ▼ 100 blades in a dispenser box, 5000 pcs/carton
- ▼ Gross Weight: 7.5KG/Carton
Net Weight: 7.0KG/Carton



■ Skin Grafts

- ▼ Surgical Skin Graft Blades are 158 mm long.
- ▼ These blades can be used with Watson, Braithwaite or Cobbett Knives.
- ▼ Available in both Stainless Steel and non-glare Carbon Steel.
- ▼ Individually wrapped in sterile aluminum foil pack by Gamma Radiation.
- ▼ 10 pcs in a dispenser box



■ **SSHIELD SAFETY SCALPEL™**
WE ONLY KEEP INNOVATING



▼ *SSHIELD Safety Scalpels™ with a sliding shield to cover the blade and easy one hand activation*

Contact Person: Mr. Shi Ye

Tel: +86 18651106599

Whatsapp: +86 18651106599

E-Mail: Szky_Ye@Hotmail.com

Website: www.medi-kyuan.com

Address: Beiqiao Town, Suzhou City, P.R China

Blood Lancet

Right blood volume less pain

Colour	Spec.	Packaging
Green	33G	
Ivory	32G	
Purple	30G	100pcs/box
Blue	28G	or 200pcs/box
Pink	26G	10000pcs/carton
Light Blue	23G	
Yellow	21G	

Contact Person: Mr.Shi Ye
Tel & Whatsapp: +86 18651106599
Email & Skype: Szky_Ye@Hotmail.com
Website: www.medi-kyuan.com
Address: Beiqiao Town, Suzhou City, P.R.China

Suzhou Kyuan Medical Apparatus Co.,LTD.

KYUAN[®]

Compatibility

fit most capillary blood lancing device



Comfort

Triple sharpened needle tips
make less pain

Safety

Sterilized by gamma radiation
One injection mold



SafetyLancet

The Only **1**-Click

Colour	Spec.	Packaging
Yellow	30G/1.8mm	
Blue	28G/1.8mm	
Purple	26G/1.8mm	50pcs/box or
Orange	23G/1.8mm	100pcs/box
Gray	21G/1.8mm	5000pcs/carton
Ivory	21G/2.2mm	
Green	18G/1.8mm	

Contact Person: Mr.Shi Ye
Tel & Whatsapp: +86 18651106599
Email & Skype: Szky_Ye@Hotmail.com
Website: www.medi-kyuan.com
Address: Beiqiao Town, Suzhou City, P.R.China

Suzhou Kyuan Medical Apparatus Co.,LTD.

KYUAN[®]

Simple

Intuitive procedure

Just twist and click

Comfort

Triple sharpened needle tips
make less pain

Ergonomic design for a
comfortable grip

Safety

Automatic retraction of lancet
avoids the possibility of injury



Ophthalmic Knife

Ophthalmic Knife

KYUAN[®]



- ◆ Cuchilletes con mango plástico desechable.
- ◆ Hay dos tipos de mango para todo tipo de cuchilletes: mango normal y corto.
- ◆ El mango normal mide 12.5cm de longitud y el mango corto mide 9cm de longitud.
- ◆ Embolsado individualmente y esterilizado por radiación gamma.
- ◆ 6 pedazos por una caja.

Slit (Phaco) Knife

Utilizada principalmente para incisiones en cirugía de facoemulsificación.



- | | | |
|---------------------|--------------------|--------------------|
| Slit Angeled 1.8mm | Slit Angeled 2.0mm | Slit Angeled 2.2mm |
| Slit Angeled 2.3mm | Slit Angeled 2.4mm | Slit Angeled 2.5mm |
| Slit Angeled 2.65mm | Slit Angeled 2.8mm | Slit Angeled 3.0mm |
| Slit Angeled 3.2mm | | |

Stab Knife

Utilizada principalmente en cirugías de catarata, glaucoma, evisceración y trasplante de córnea.



- | | | | |
|----------|------------|----------|----------|
| Stab 15° | Stab 22.5° | Stab 30° | Stab 45° |
|----------|------------|----------|----------|

Implant Knife

Utilizada principalmente en cirugías de catarata y trasplante de córnea.



- | | |
|-------------------------------------|-------------------------------------|
| Angled 3.5mm (For Fordable IOL) | Angled 4.1mm (For Fordable IOL) |
| Angled 5.2mm (For Fordable IOL) | Angled 5.5mm (For Non-Fordable IOL) |
| Angled 6.0mm (For Non-Fordable IOL) | |

Crescent Knife

Utilizada principalmente para incisiones en cirugía de Trabeculectomia y Catarata (EECC).



Angled 2.0mm

Angled 2.2mm

MVR

Utilizada principalmente para incisiones en cirugía de esclerotomia.



- | | | |
|--------------|-------------|-------------|
| 19G (1.00mm) | 20G(0.90mm) | 21G(0.80mm) |
| 22G(0.70mm) | 23G(0.60mm) | |

KYUAN[®]

SSHIELD SAFETY OPHTHALMIC KNIFE[™]
WE ONLY KEEP INNOVATING



- ▼ *SSHIELD Safety Ophthalmic Knife[™] with a sliding shield to cover the blade and easy one hand activation*

Contact Person: Mr. Shi Ye

Tel: +86 18651106599

Whatsapp: +86 18651106599

E-mail and Skype: Szky_Ye@Hotmail.com

Website: www.medi-kyuan.com

Address: Beiqiao Town, Suzhou City, P.R China